# CombiScreen® mALB / CREA

REF	Product	Content
94025	CombiScreen® mALB / CREA	25 Urine test strips

## INTENDED USE

The CombiScreen® mALB / CREA urine test strip is intended to be used for the semi-quantitative determination of albumin and creatinine in human urine<sup>1</sup>. This product can be used for screening samples for microalbuminuria<sup>2</sup>. The mALB / CREA test strip is designed as IVD and only for professional use. It should only be done by medically trained personnel. Evaluation is performed visually only.

## SUMMARY AND EXPLANATION

Microalbuminuria is a symptom of beginning nephropathy and can be a signal for cardio-vascular diseases<sup>3, 4</sup>. The albumin concentration in the urine varies depending on the voided urine volume<sup>5</sup>. Creatinine excretion into the urine is usually constant. Using the albumin to creatinine ratio it is possible to correct the volume effect and to use spontaneous voided instead of 24-hour urine for the diagnosis of microalbuminuria<sup>6</sup>.

## TEST PRINCIPLE

Albumin: The test is based on the "protein error" principle of a tetrabromophenol sulfone phthalein derivative as indicator<sup>7</sup>. Under acidic conditions, binding of the dye to albumin leads to a color change from light to dark turquoise. Creatinine: The test is based on the peroxidase-like activity of a copper-creatinine-complex. This complex catalyzes the color reaction from light-green to dark blue-green<sup>8</sup>.

## REAGENTS

Albumin test pad: Tetrabromophenol sulfone phthalein derivative 1,6%  
Creatinin test pad: Copper sulfate 1,5%; Cumolhydroperoxide 4%; Tetramethylbenzidine 1,7%

## Warnings and precautions

Only for in-vitro diagnostic use.  
The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. In case any serious incident has occurred in relation to the device, please report to the manufacturer and, if applicable, the competent authority of the country in which the user and/or the patient is established.  
Do not touch the test pads! Avoid ingestion, eye contact and contact with mucous membranes. Keep away from children.  
For safe handling of urine test strips and for avoiding contact with potentially infectious substances, please follow the general working instructions for laboratories.  
Dispose of the used test strips in accordance with the valid safety regulations.  
External influences such as humidity, light and temperatures can cause a discoloration of reagent areas and may indicate deterioration. Do not use discolored urine test strips.

## Storage and stability

Keep diagnostic test strips protected from direct sunlight and humidity. Store the tubes in a cool and dry place (storage temperature 2–22 °C). The test strips can be used until the given expiry date if stored under appropriate conditions.

## COLLECTION AND PREPARATION

It is recommended to use only fresh, native, well-mixed and non-centrifuged urine. Test within 2 hours, protect samples from light. Morning urine is preferable. Collection of 24-hour urine is not required<sup>9</sup>. If direct testing is not applicable, store sample at 2–4 °C. Warm up to room temperature (15–25 °C) before testing. Collection tubes must be clean, dry and free from detergents, biocides or disinfectants. Do not add preservatives.

## PROCEDURE

- Remove only the necessary number of test strips for measurement immediately before use and directly close the tube after removal tightly with the stopper.
- Dip the test strip shortly (approx. 1–2 seconds) into the well-mixed urine. Make sure that both test pads are immersed in the sample. Remove excess urine from the strip by wiping the edge of the strip on the urine container or on absorbent paper.
- To prevent interaction from adjacent test areas, hold the strip in a horizontal position during incubation.
- Compare the reagent areas on the strip with the corresponding color chart on the container 60 seconds after immersion. Visual evaluation should be carried out in diffuse daylight. Coloration only on the rim of the test pad or after more than 2 minutes after immersion is without meaning and should not be used for interpretation.
- Assess the albumin/creatinine ratio using Table 1.

## Materials provided

CombiScreen® mALB / CREA Test Strips  
Materials required but not provided  
Clean and dry collection tubes, free from detergents, biocides or disinfectants.

## QUALITY CONTROL

Test strip performance can be checked with the CombiScreen® controls or other appropriate control solutions. The use of quality control solutions should follow the internal guidelines of the laboratory and the local regulations. It is recommended to perform control measurements after opening a new test strip tube or use of a new test strip batch. Each laboratory is recommended to establish its own quality control standards.

## RESULTS

After identification of the albumin and creatinine concentration the albumin/creatinine ratio can be assessed as "normal", "abnormal" or "high abnormal" according to the interpretation table on the tube label.  
[X]: Sample is too diluted. Repeat the test with a new sample.  
The ratio can also be given in milligram albumin per gram creatinine (mg/g) or as milligram albumin per millimol creatinine (mg/mmol).

**Table 1: Interpretation table for the albumin/creatinine ratio**

Albumin	Creatinine				
	A (10 mg/dL) (0,9 mmol/L)	B (50 mg/dL) (4,4 mmol/L)	C (100 mg/dL) (8,8 mmol/L)	D (200 mg/dL) (17,7 mmol/L)	E (300 mg/dL) (26,5 mmol/L)
1 (10 mg/L)	X	Normal	Normal	Normal	Normal
2 (30 mg/L)	High Abnormal	Abnormal	Normal	Normal	Normal
3 (80 mg/L)	High Abnormal	Abnormal	Abnormal	Abnormal	Normal
4 (150 mg/L)	High Abnormal	High Abnormal	Abnormal	Abnormal	Abnormal
5 (500 mg/L)	High Abnormal	High Abnormal	High Abnormal	Abnormal	Abnormal

## Example for result interpretation

If a urine sample contains 10 mg/L albumin and 100 mg/dL creatinine, the albumin/creatinine ratio is classified as normal. A test result of 150 mg/L albumin and 200 mg/dL creatinine is classified as abnormal (microalbuminuria).

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## LIMITATIONS OF THE PROCEDURE

Detergents, cleaning agents, disinfectants and preservatives may lead to false values for the albumin or creatinine concentration. Different urine contents, especially high concentrations of hemoglobin, riboflavin or bilirubin, can lead to atypical coloration on the test pads.

In order to establish a final diagnosis and prescribe an appropriate therapy, the results obtained with test strips should be evaluated in combination with other medical results and the patient's medical history.  
The effect of medicaments, drugs or their metabolic products on the test is not known in all cases. In case of doubt it is therefore recommended to repeat the test after discontinuation of the medication. However, stopping medication should only be done after respective instruction of the doctor.  
Due to the fact that the content of the urine is not constant (e.g. content of activators or inhibitors which may vary from sample to sample, changing ion concentration), the conditions of the reaction are not always the same which may lead to variations of the intensity and the color in rare cases.

## EXPECTED VALUES

**Table 2: Expected values of the albumin/creatinine ratio:**

Classification	Conv. unit (mg/g)	SI unit (mg/mmol)	Indication for:
Normal	≤ 30	≤ 3.4	
Abnormal	31–299	3.5–33.8	Microalbuminuria
High abnormal	≥ 300	≥ 33.9	Macroalbuminuria, Proteinuria

## PERFORMANCE CHARACTERISTICS

The performance characteristics of the CombiScreen® mALB / CREA urine test strips have been determined on the basis of analytical performance studies. Standard routine urine samples from doctors' offices were used for performance studies. These samples were analyzed with CombiScreen® mALB / CREA urine test strip and with an immunological albumin method and an enzymatic creatinine assay<sup>10</sup>. Accuracy of the urine test strips was characterized by its agreement with the comparative methods. Table 3 shows the results of the performance studies.

**Table 3: Performance of the CombiScreen® mALB / CREA urine test strip relative to comparative methods.**

Result	a) Accuracy	b) positive percent agreement (PPA) / diagnostic sensitivity	c) negative percent agreement / diagnostic specificity
Albumin/Creatinine Ratio	85.7 % (n = 280)	93.1 % (n = 280)	91.1 % (n = 280)

## Analytical sensitivity:

The test reports values from minimum 10 mg/L albumin and 10 mg/dL creatinine.

## Precision:

20 parallel determinations with three different lots of CombiScreen® mALB / CREA urine test strips were performed. Each lot was tested on two urine specimen (Within-Run). For reproducibility data (Day-to-Day), single detections were performed on 20 days with three lots of CombiScreen® mALB / CREA urine test strips. The precision was calculated from the number of matched readings on the same specimen.

Albumin concentration	Within-Run	Day-to-Day	Creatinine concentration	Within-Run	Day-to-Day
30 mg/L	100 %	100 %	50 mg/dL	100 %	100 %
150 mg/L	100 %	95 %	200 mg/dL	100 %	100 %

## Analytical specificity/Interferences

The following substances do not affect the test system up to the specified concentration:  
Glucose 1500 mg/dL; Bilirubin 1 mg/dL; Ascorbic acid 80 mg/dL; Urobilinogen 12 mg/dL; Sodium nitrite 1 mg/dL; Leukocytes 500 mg/dL; Hemoglobin 0.87 mg/dL; Lithium acetoacetate 300 mg/dL; Acetylsalicylic acid 60 mg/dL; Caffeine 60 mg/dL; Riboflavin 10 mg/dL; Creatine 10 mg/dL; Formaldehyde 0.1%; Fructose 100 mg/dL; Galactose 80 mg/dL; Uric acid 50 mg/dL; Urea 400 mg/dL; Lactose 10 mg/dL; Human IgG 5 mg/dL; Myoglobin 0.26 mg/dL. Alkaline samples can produce false-positive results on the albumin test panel.

## Measuring range:

The color changes of the test pads correspond to the following concentrations:  
Albumin: 10, 30, 80, 150, 500 mg/L.  
Creatinine: 10 (0,9), 50 (4,4), 100 (8,8), 200 (17,7), 300 mg/dL (26,5 mmol/L).

## LITERATURE

- Parsons M, Newman DJ, Puglia M, Newall RG, and Price CP; Performance of a reagent strip device for quantitation of the urine albumin: creatinine ratio in a point of care setting. Clin Nephrol. 1999 Apr;51(4):220-7.
- Sasaki M, Puglia MJ, Parker DR, Kuromoto K, Furukawa I, and Konishi I. Measurement of the albumin content of urinary protein using dipsticks. J Clin Lab Anal. 1999;13(5):246-50.
- Kaplan NM. Microalbuminuria: a risk factor for vascular and renal complications of hypertension. Am J Med. 1992 Apr 27;92(4B):8S-12S.
- Apostolovic S, Stanojevic D, Djordjevic V, Jankovic Tomasevic, R, Salinger Martinovic S, Djordjevic Radojkovic D, Koracevic G, Stojanovic I, and Pavlovic M. Prognostic significance of microalbuminuria in patients with acute myocardial infarction. Clin. Lab. 57; 229-235; 2011
- Puglia MJ, Lott JA, Clark LW, Parker DR, Wallace JF, Willis TW. Comparison of urine dipsticks with quantitative methods for microalbuminuria. Eur J Clin Chem Clin Biochem. 1997 Sep;35(9):693-700.
- Position Statement: Nephropathy in Diabetes. Diabetes Care, 27, S79-S83, 2004
- Puglia MJ, Lott JA, Proffitt JA, and Cast TK. High-sensitivity dye binding assay for albumin in urine. J Clin Lab Anal. 1999;13(4):180-7.
- Puglia MJ, Lott JA, Wallace JF, Cast TK, Bierbaum LD. Assay of creatinine using the peroxidase activity of copper-creatinine complexes. Clin Biochem. 2000 Feb;33(1):63-70.
- Guy M, Borzomato JK, Newall RG, Kalra PA, and Price CP . Protein and albumin-to-creatinine ratios in random urines accurately predict 24 h protein and albumin loss in patients with kidney disease. Ann Clin Biochem. 2009;46: 468-76.
- Burtis CA, Ashwood ER and Bruns DE (eds): Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (5th edition), Elsevier, St. Louis, USA, 2012, 2238 pp, 909 illustrations. ISBN: 978-1-4160-6164-9

## Symbols

- = In vitro diagnostics product
- = The product complies with European legislation
- = Follow the instructions for use!
- = Use by
- = Permitted storage temperature range
- = Only single use
- = Batch identification number
- = Item number
- = Manufacturer

Analyticon Biotechnologies AG  
Am Muehlenberg 10  
35104 Lichtenfels  
Germany  
[www.analyticon-diagnostics.com](http://www.analyticon-diagnostics.com)

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## Literatur

- Parsons M, Newman DJ, Puglia M, Newall RG, and Price CP; Performance of a reagent strip device for quantitation of the urine albumin: creatinine ratio in a point of care setting. Clin Nephrol. 1999 Apr;51(4):220-7.
- Sasaki M, Puglia MJ, Parker DR, Kuromoto K, Furukawa I, and Konishi I. Measurement of the albumin content of urinary protein using dipsticks. J Clin Lab Anal. 1999;13(5):246-50.
- Kaplan NM. Microalbuminuria: a risk factor for vascular and renal complications of hypertension. Am J Med. 1992 Apr 27;92(4B):8S-12S.
- Apostolovic S, Stanojevic D, Djordjevic V, Jankovic Tomasevic, R, Salinger Martinovic S, Djordjevic Radojkovic D, Koracevic G, Stojanovic I, and Pavlovic M. Prognostic significance of microalbuminuria in patients with acute myocardial infarction. Clin. Lab. 57; 229-235; 2011
- Puglia MJ, Lott JA, Clark LW, Parker DR, Wallace JF, Willis TW. Comparison of urine dipsticks with quantitative methods for microalbuminuria. Eur J Clin Chem Clin Biochem. 1997 Sep;35(9):693-700.</

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- Parsons M, Newman DJ, Pugia M, Newall RG, and Price CP: Performance of a reagent strip device for quantitation of the urine albumin: creatinine ratio in a point of care setting. Clin Nephrol. 1999 Apr;51(4):220-7.
- Sasaki M, Pugia MJ, Parker DR, Kuromoto K, Furukawa I, and Konishi I. Measurement of the albumin content of urinary protein using dipsticks. J Clin Lab Anal. 1999;13(5):246-50.
- Kaplan NM. Microalbuminuria: a risk factor for vascular and renal complications of hypertension. Am J Med. 1992 Apr 27;92(4B):85-125.
- Apostolovic S, Stanojevic D, Djordjevic V, Jankovic Tomasevic, R, Salinger Martinovic S, Djordjevic Radojkovic D, Koracevic G, Stojanovic I, and Pavlovic M. Prognostic significance of microalbuminuria in patients with acute myocardial infarction. Clin. Lab. 57; 229-235; 2011
- Pugia MJ, Lott JA, Clark LW, Parker DR, Wallace JF, Willis TW. Comparison of urine dipsticks with quantitative methods for microalbuminuria. Eur J Clin Chem Clin Biochem. 1997 Sep;35(9):693-700.
- Position Statement: Nephropathy in Diabetes. Diabetes Care, 27, S79-S83, 2004
- Pugia MJ, Lott JA, Proffitt JA, and Cast TK. High-sensitivity dye binding assay for albumin in urine. J Clin Lab Anal. 1999;13(4):180-7.
- Pugia MJ, Lott JA, Wallace JF, Cast TK, Bierbaum LD. Assay of creatinine using the peroxidase activity of copper-creatinine complexes. Clin Biochem. 2000 Feb;33(1):63-70.
- Guy M, Borzomato JK, Newall RG, Kalra PA, and Price CP . Protein and albumin-to-creatinine ratios in random urines accurately predict 24 hr protein and albumin loss in patients with kidney disease. Ann Clin Biochem. 2009;46: 468-76.
- Burtis CA, Ashwood ER and Bruns DE (eds): Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (5th edition), Elsevier, St. Louis, USA, 2012, 2238 pp, 909 illustrations. ISBN: 978-1-4160-6164-9

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